






HIGH RESOLUTION MANOMETRY CATHETERS

Order No.	Description	Price
<p>HRM36</p>	 <p>Solid State 36 Channel Catheter, circumferential sensors Diameter 12/16 french Can be used with ISOLAB HR or MALT</p>	
<p>HRM36DIR</p>	 <p>Solid State 36 Channel Catheter, directional sensors Diameter 8 french Can be used with ISOLAB HR or MALT</p>	
<p>HRM36IMP16</p>	 <p>Solid State 36 Channel Catheter, circumferential sensors 16 Impedance sensors (15 channels) Diameter 12/16 french Can be used with ISOLAB HR or MALT</p>	