






Order No.	Description	Price
HRM36IMP16DIR	 <p>Solid State 36 Channel Catheter, directional sensors 16 Impedance sensors (15 channels) Diameter 8 french Can be used with ISOLAB HR or MALT</p>	
HRM8	 <p>Solid state 8 Channel Catheter, circumferential sensors, with integrated balloon entry Diameter 12/16 french Can be used with ISOLAB HR or MALT</p>	
HRMAFT3D	 <p>Solid State 33 Channel Catheter, directional sensors, with integrated balloon entry Diameter 12 french Can be used with ISOLAB HR or MALT</p>	

Notice:

In order to use the Chicago Classification for esophageal manometry the HRM catheter to be used should have circumferential sensors.