






DISPOSABLE IMPEDANCE pH-PROBES WITH INTERNAL REFERENCE, ANTIMONY, NON-INFUSED

Order No.	Description	
<p>SI-ZNIS-L</p>	 <p> -15 cm -10 cm -5 cm 0 cm 5 cm 10 cm 15 cm 20 cm SI-ZNIS-L ■ impedance ● pH ≥18 cm oesophageal length 6-Channel Impedance / 1-Channel Antimony pH Batteries Included </p>	
<p>SI-ZNIS-M</p>	 <p> -15 cm -10 cm -5 cm 0 cm 5 cm 10 cm 15 cm 20 cm SI-ZNIS-M ■ impedance ● pH 15-18 cm oesophageal length 6-Channel Impedance / 1-Channel Antimony pH Batteries Included </p>	
<p>SI-ZNIS-S</p>	 <p> -15 cm -10 cm -5 cm 0 cm 5 cm 10 cm 15 cm 20 cm SI-ZNIS-S ■ impedance ● pH <15 cm oesophageal length 6-Channel Impedance / 1-Channel Antimony pH Batteries Included </p>	