


Accessories Endo- Urology

Pelvitec



Pelvitec

Rotterdamseweg 183-C
2629 HD Delft
The Netherlands


Phone: 0031 (0)10 3032290

e-mail: info@pelvitec.nl

Website: www.pelvitec.nl

Urological guidewires


Improve safety, Smooth, Easy navigation




Hydrophilic coating on guidewire significantly reduces surface friction.

Hybrid ball tip guidewire

PTFE ball tip guidewire




Ball and atraumatic tip is designed to reduce tip impactions



Our guidewires

Description	Shape tip	Outer ϕ (inch)	Length (mm)
Hybrid guidewire	Straight	0.035"	1500
Hybrid ball tip guidewire			
PU guidewire			
PTFE guidewire	Shape	0.035"	1500
PTFE ball tip guidewire			

Ureteral Access Sheaths



Sizes Access sheaths

Sheath Inner Diameter(Fr)	Sheath Length(cm)
10	35
10	45
12	35
12	45
14	35
14	45

Urological Stone Extraction Baskets

Various types

3-wire and 4-wire design

3-wire design

3-wire basket designed to capture small stones.

4-wire design


Tipless design

Tipless design for easy access and capture stone in renal calys

Models and sizes stone baskets

Basket diameter (mm)	Article No.	Working length (mm)	Catheter OD (Fr)
8	4-wire design	1300	1.9
11			
14			
17	3-wire design		
20			

Ureteral Dilation Balloon Catheters



Inflator



Sizes dilation balloon catheters


Ureteroscope type (M type)

Catheter OD (FR)	Balloon diameter (mm)	Balloon length (mm)	Catheter length (mm)
2.4	4 5 7	60	1500

Wireguided type (N type)

Catheter OD (FR)	Balloon diameter (mm)	Balloon length (mm)	Catheter length (mm)
6	8 10	60	1500

Re-usable Biopsy Forceps, 120 cm, 3 FR




Extremely sharp and specially designed cups to realize a cutting and not a crushing biopsy

Outstanding flexibility and operational reliability



360° axial rotation


Re-usable Forceps



Alligator forceps, 1,0 mm 3FR, 120 cm










Rat-toothed forceps, 1,0 mm 3FR, 120 cm



Long-oval cups serrated forceps, 1,0 mm 3FR, 120 cm

Order information

Article number	Description
	Guidewires
18-UG3515USO	Guide Wire, PU, Nitinol, 0.035i x 150 cm
18-UG3515BSO	Guide Wire, PU+PTFE hybrid, Nitinol, Hydrophilic Tip, 0.035i x 150 cm
	Guidewires ball tip
18-UG3515TSO	Guide Wire, Ball Tip, PTFE, Nitinol, 0.035i x 150 cm
18-UG3815CSO	Guide Wire, Ball Tip, PU+PTFE hybrid, Hydrophilic Tip, 0.035i x 150 cm
	Access sheaths
18-UA1035	Access sheath, hydrophilic, 10/12 F 35cm
18-UA1045	Access sheath, hydrophilic, 10/12 F 45cm
18-UA1135	Access sheath, hydrophilic, 11/13 F 35cm
18-UA1145	Access sheath, hydrophilic, 11/13 F 45cm
18-UA1235	Access sheath, hydrophilic, 12/14 F 35cm
18-UA1245	Access sheath, hydrophilic, 12/14 F 45cm
18-UA1335	Access sheath, hydrophilic, 13/15 F 35cm
18-UA1345	Access sheath, hydrophilic, 13/15 F 45cm
18-UA1435	Access sheath, hydrophilic, 14/16 F 35cm
18-UA1445	Access sheath, hydrophilic, 14/16 F 45cm
	Stone Baskets van het type 3-wire open
18-UST0810	Stone Basket, 3-Wire, 1,9 F, 8 mm x 130 mm
18-UST1110	Stone Basket, 3-Wire, 1,9 F, 11 mm x 130 cm
18-UST1410	Stone Basket, 3-Wire, 1,9 F, 14 mm x 130 cm
18-UST1710	Stone Basket, 3-Wire, 1,9 F, 17 mm x 130 mm
18-UST2010	Stone Basket, 3-Wire, 1,9 F, 20 mm x 130 cm

Article number	Description
	Stone Baskets van het type Dormia
18-USF0810	Stone Basket, Tipless, 4-Wire, 1,9 F, 08 mm x 130cm
18-USF1110	Stone Basket, Tipless, 4-Wire, 1,9 F, 11 mm x 130 cm
18-USF1410	Stone Basket, Tipless, 4-Wire, 1,9 F, 14 mm x 130 cm
18-USF1710	Stone Basket, Tipless, 4-Wire, 1,9 F, 17 mm x 130 cm
18-USF2010	Stone Basket, Tipless, 4-Wire, 1,9 F, 20 mm x 130 cm
	Ureteral Dilation Balloon Catheter
18-UBUN08060	Dilation Cath. 7F, for Guide Wire, 8x60 mm balloon
18-UBUN10060	Dilation Cath. 7F, For Guide Wire, 10x60 mm balloon
	Ureteral Dilation Balloon Catheter for rigid scope
18-UBMN04060	Dilation Cath. 2.4F, Rigid Scope, 4x60 mm balloon
18-UBMN05060	Dilation Cath. 2.4F, Rigid Scope, 5x60 mm balloon
18-UBMN07060	Dilation Cath. 2.4F, Rigid Scope, 7x60 mm balloon
	Re-usable Forceps
18-6100037	Biopsy forceps,flexible,3Fr/CH,narrowing distal shaft/spring
18-6100032	Alligator forceps, 1.0mm 3Fr/CH, narrow.distal shaft/spring
18-6100034	Rat-toothed forceps,1.0mm 3Fr/CH, narrow.distal shaft/spring
18-6100036	Long-oval cups serrated-forceps,1.0mm 3Fr/CH/spring

Pelvitec